

## **BIOPSY PADS**



CODE	COLOUR	PACKAGING
07-7290	Blue	500 pz
07-7290/5	Blue	5000 pz
07-7290/15	Blue	15000 pz
07-7291	Black	500 pz



In Vitro Diagnostic – Medical Device



Manufacturer: Bio-Optica Milano S.p.A.



Disposable Product

Date of emission: 24/03/2020

Rev. 03

Special pads for the embedding of biopsies and small anatomical specimens, to hold them in place and prevent lost during processing.

**TECHNICAL SPECIFICATIONS**

Polyurethane foam guarantees high solvents resistance and its structure with open cells allows fluids (alcohols, solvents and paraffin) to transfer safely through the tissue during the processing.

Filters can be autoclaved (max. Temperature: 134°C). For sterilization, it is enough a temperature of 121°C.

Filters are packed in hermetic plastic bags, to maintain them clean and intact until the use.

Technical data	Chemical composition	Polyester - polyurethane foam (moltoprene)
	Dimensions	2.5 x 3 cm (thikness: 2 mm)
Packaging	Primary container	Plastic bag
	Secondary container	Carton box
Storage	Storing conditions	None
	Validity	Not applicable to this kind of product
Warning and precautions	Product classification	For specialized technical operators. Product not classified as dangerous.
	Disposal	Dispose in conformity with local laws.

Date of emission: 24/03/2020

Rev. 03